BG051PO001-3.3.06 -0059

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ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

Quick-RNA™ MiniPrep

Catalog Nos. **R1054 & R1055**

Highlights

- High-quality total RNA from a wide range of samples.
- You can opt to isolate small and large RNAs in separate fractions.
- *DNA-free* RNA is ready for use in any downstream application. *DNase I included.*

Contents

Product Contents	1
Specifications	1
Product Description	2
Buffer Preparation	3
Protocols	3, 4
Appendices	5
Ordering Information	6
Related Products	7-8

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Satisfaction of all Zymo Research products is guaranteed. If you are dissatisfied with this product please contact us.

Some difficult-to-lyse samples may require mechanical or enzymatic homogenization. For assistance, contact us at tech@zymoresearch.com.

Use the **Quick-RNA™ MicroPrep** (Cat. Nos. R1050, R1051) for up to 10 µg RNA from 1-10⁸ cells.

Product Contents

Quick-RNA™ MiniPrep (Kit Size)	R1054 (50 Preps.)	R1055 (200 Preps.)
RNA Lysis Buffer	50 ml	2x 100 ml
RNA Prep Buffer	25 ml	100 ml
RNA Wash Buffer¹ (concentrate)	24 ml	2x 48 ml
DNase/RNase-Free Water	6 ml	3x 10 ml
DNase I Set² DNase I (250 U) & 10x DNase I Reaction Buffer (1 ml)	1 set	4 sets
Spin-Away™ Filters	50	200
Zymo-Spin™ IIICG Columns	50	200
Collection Tubes	100	400
Instruction Manual	1	1

Note – Integrity of kit components is guaranteed for up to one year from date of purchase. Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

Storage Temperature - Store all kit components (*i.e.*, buffers, columns) at room temperature. Store reconstituted DNase I at -20 °C.

¹ Add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate or 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate before use.

² Add 275 µl **DNase/RNase-Free Water** per vial to reconstitute the lyophilized **DNase I** (E1009) at 1 U/µl. Mix by gentle inversion. Store aliquots at -20 °C.

Specifications

- **Sample Sources** – Cells or tissue samples, yeast, plant or bacteria. Compatible with DNA/RNA Shield™ (Cat. No. R1100) and RNAlater™.
- **Sample Storage** – Samples homogenized in RNA Lysis Buffer are stable and can be stored frozen prior to purification.
- **Sample Size** – Up to 10⁷ cells or 50 mg tissue.
- **RNA Purity** – High quality RNA ($A_{260}/A_{280} > 1.8$, $A_{260}/A_{230} > 1.8$) suitable for all downstream RNA-based manipulations.
- **RNA Recovery** – Up to 100 µg RNA can be eluted into ≥30 µl RNase-free water allowing for a highly concentrated sample.
- **RNA Storage** – RNA is eluted with RNase-free water and can be stored frozen. RNase inhibitors can be included for prolonged storage.
- **Equipment Needed** – Microcentrifuge.

Note - ™ Trademarks of Zymo Research Corporation. This product is for research use only and should only be used by trained professionals. It is not intended for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility. RNAlater™ is a trademark of Ambion, Inc.

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

Product Description

The **Quick-RNA™ MiniPrep** kit is an innovative product designed for the easy, reliable, and rapid isolation of DNA-free RNA from a wide range of cell (*up to 10⁷*) and tissue samples (*up to 50 mg*). The procedure combines a unique buffer system with *Clean-Spin* column technology to yield high quality total RNA (*including small RNAs 17-200 nt*) in about 10 minutes.



The procedure is simple. Add the provided **RNA Lysis Buffer** to a sample, and then purify the RNA using the **Zymo-Spin™ Columns**. The result is highly-concentrated, *DNA-free* RNA that is suitable for RT-PCR, hybridization, sequencing *etc.* In addition, the kit can be used for the enrichment of small and large RNAs into separate fractions (Appendix A, page 5).

For Assistance, please contact Zymo Research Technical Support at 1-888-882-9682 or e-mail tech@zymoresearch.com.

Notes:

Use the **Direct-zol™ RNA MiniPrep** (Cat. Nos. R2050, R2051, R2052, R2053) for isolation of RNA directly (without phase separation) from samples in Trizol®, *etc.*

Use the **DNA/RNA Shield™** (Cat. Nos. R1100-50, R1100-250) for safe sample storage and transport at ambient temperatures.

The **Quick-RNA™ MiniPrep** yields high quality total RNA. High levels of genomic DNA contamination are present in the preps from Suppliers Q & P but not with the **Quick-RNA™ MiniPrep**. Total RNA was isolated from human epithelial cells (sans DNase treatment).

RNA isolated with the **Quick-RNA™ MiniPrep** is DNA-free. Samples isolated with Supplier Q's kit are provided for comparison. Total RNA was isolated from 10⁶ human epithelial cells (with in-column DNase treatments for both kits). Each amplification curve represents an average of three independent isolation experiments.

The **Quick-RNA™** kits yield high quality RNA as indicated by the RIN (RNA Integrity Number; 2200 TapeStation, Agilent).

RNA MiniPrep Kit Comparison

	Quick-RNA™	Supplier Q
Small RNA (≥17 nt) recovery	YES	NO
DNase I included	YES	NO
gDNA removal column included	YES	NO

ZYMO RESEARCH CORP.

Ensure the RNA isolation procedure is performed in an RNase-free environment.

Buffer Preparation

- ✓ Before starting, add 96 ml 100% ethanol (104 ml 95% ethanol) to the 24 ml **RNA Wash Buffer** concentrate (R1054) or 192 ml 100% ethanol (208 ml 95% ethanol) to the 48 ml **RNA Wash Buffer** concentrate (R1055).
- ✓ Add 275 µl **DNase/RNase-Free Water** per vial to reconstitute the lyophilized **DNase I** (E1009) at 1 U/µl. Mix by gentle inversion. Store aliquots at -20 °C.

Protocols

The RNA isolation consists of three steps: (I) *Sample Lysis/Homogenization*, (II) *Sample Clearing and gDNA Removal* and (III) *RNA Purification*.

All steps should be performed at room temperature (20-30 °C).

I. Sample Lysis/Homogenization

Recommended RNA Lysis Buffer volumes

RNA Lysis Buffer	300 µl	600 µl
Cells	Up to 5 x 10 ⁶	>5 x 10 ⁶
Tissue	<20 mg	≤50 mg

Adherent Cells

Lyse cells directly in the culture container by removing liquid medium and adding **RNA Lysis Buffer** directly to the monolayer.

Cells in Suspension

Pellet cells (≤500 x g), remove the supernatant completely then resuspend the cell pellet in **RNA Lysis Buffer**. Vortex briefly.

Tissue and Tough-to-Lyse Samples

Fresh or frozen tissue (animal, plant, insect, yeast or bacteria) can be mechanically homogenized (e.g., **ZR BashingBead™ Lysis Tubes**) directly in the **RNA Lysis Buffer**.

Alternatively, tough-to-lyse tissue samples can be Proteinase K treated (Appendix B, page 5).

Liquids/Reaction Clean-up

DNase-treated RNA, labeling and *in vitro* transcription reactions can be processed directly by adding 4 volumes of **RNA Lysis Buffer** to each volume of sample (4:1) then mixing well.

Samples in DNA/RNA Shield™

Bring samples homogenized and stored in **DNA/RNA Shield™** to room temperature (20-30 °C). Then add 1 volume **RNA Lysis Buffer** (1:1), mix and proceed with Sample Clearing step.

Samples in DNA/RNA Shield™ can be Proteinase K treated (Appendix B, page 5).

Samples in RNA/later™

To process cells or liquids in **RNA/later™** (without reagent removal): Add 1 volume of RNase-free water or PBS to the sample (1:1). Then add 4 volumes **RNA Lysis Buffer** (4:1) and mix.

Alternatively, remove the **RNA/later™**, then proceed with Sample Lysis/Homogenization according to the sample type.

Notes:

Samples homogenized in **RNA Lysis Buffer** can be stored frozen for processing at a later time.

ZR Bashing Bead™ Lysis Tubes are available separately (Cat. Nos. S6002, S6003).

Processing plant tissue and other samples containing polyphenolics, humic acids, melanin, etc. may require use of the **OneStep™ PCR Inhibitor Removal Kit** (Cat. No. D6030).

Use the **DNA/RNA Shield™** (Cat. Nos. R1100-50, R1100-250) for safe sample storage and transport at ambient temperatures.

ZYMO RESEARCH CORP.

Appendix A: Purification of Small and Large RNAs into Separate Fractions

This procedure is compatible with animal cell inputs (up to 10⁶) or previously isolated RNA only.

All centrifugation steps should be performed between 10,000-16,000 x g. This protocol requires two columns (per prep).

- Mix an equal volume of **RNA Lysis Buffer** and ethanol (95-100%).
Example: Mix 50 µl buffer and 50 µl ethanol.
- Add 2 volumes of the buffer/ethanol to an RNA sample¹ or 300 µl buffer/ethanol to a cell pellet and mix.
Example: Mix 100 µl buffer/ethanol and 50 µl sample.
- Transfer the mixture² to the **Zymo-Spin™ Column** and centrifuge for 30 seconds. **Save the flow-through!**

Column: RNAs >200 nt


Flow-through: RNAs 17-200 nt

- Continue to step 5.
- Add 400 µl **RNA Prep Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
- Add 700 µl **RNA Wash Buffer** to the column and centrifuge for 30 seconds. Discard the flow-through.
- Add 400 µl **RNA Wash Buffer** and centrifuge the column for 2 minutes to ensure complete removal of the wash buffer.
- Place the column into an RNase-free tube (not provided). Add ≥30 µl **DNase/RNase-Free Water**³ directly to the column matrix, then centrifuge at top speed for 30 seconds.
Eluted RNA can be used immediately or stored frozen.

Add 1 volume ethanol and mix.

Example: Add 150 µl ethanol to 150 µl flow-through.

Transfer the mixture to a new column and centrifuge for 30 seconds. Discard the flow-through.

Total RNA (>17 nt), large (>200 nt) or small RNAs (17-200 nt) are effectively partitioned and purified with the **Quick-RNA™** kit.

Notes:

¹ Adjust the sample volume to 50 µl (minimum).

² Zymo-Spin™ columns may be reloaded to process samples >700 µl.

³ 2X Digestion Buffer (Cat. No. D3050-1-5 and D3050-1-20).

⁴ Proteinase K (Cat. No. D3001-2-5 and D3001-2-20).

One unit of enzyme will hydrolyze urea-denatured hemoglobin to produce 1.0 µmole of tyrosine per minute at pH 7.5 at 37°C.

Appendix B: Proteinase K Digestion

Example: up to 5 mg solid tissue or 10⁶ animal cells in DNA/RNA Shield™
2X Digestion Buffer³
Proteinase K⁴

95 µl
95 µl
≥6 U

Prepare a Proteinase K reaction mix (see example above, scale-up as necessary). Incubate at 55°C for 30 minutes (e.g., pelleted white blood cells) or 1-3 hours (solid tissue). Then add 1 volume **RNA Lysis Buffer** and proceed to Sample Clearing and gDNA Removal (page 4).

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 • Toll Free: (888) 882-9682 • Fax: (949) 266-9452 • info@zymoresearch.com • www.zymoresearch.com

Ordering Information

Product Description	Input	Binding	Catalog No.	Kit Size
Quick-RNA™ MicroPrep	~1-10 ⁶ cells	~10 µg	R1050	50 Preps.
			R1051	200 Preps.
Quick-RNA™ MiniPrep	~10 ² -10 ⁷ cells	~100 µg	R1054	50 Preps.
			R1055	200 Preps.
Quick-RNA™ MidiPrep	~10 ⁶ -10 ⁸ cells	~1 mg	R1056	25 Preps.
ZR-96 Quick-RNA™	~1-10 ⁶ cells	~10 µg/well	R1052	2x 96 Preps.
			R1053	4x 96 Preps.

For Individual Sale	Catalog No.	Amount
RNA Lysis Buffer	R1060-1-50	50 ml
	R1060-1-100	100 ml
RNA Prep Buffer	R1060-2-10	10 ml
	R1060-2-25	25 ml
	R1060-2-100	100 ml
	R1003-3-6	6 ml
RNA Wash Buffer (concentrate)	R1003-3-12	12 ml
	R1003-3-24	24 ml
	R1003-3-48	48 ml
	DNase I (lyophilized) (250 U supplied with 10x DNase I Reaction Buffer)	E1009
Spin-Away™ Filter	C1006-50-F	50
	C1006-250-F	250
Zymo-Spin™ IICG Column	C1006-50-G	50
	C1006-250-G	250
Collection Tube	C1001-50	50
	C1001-500	500
	C1001-1000	1000
DNase/RNase-Free Water	W1001-1	1 ml
	W1001-6	6 ml
	W1001-10	10 ml

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

DNA PURIFICATION

What is Clean-Spin™ Technology?
 The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Purify DNA from PCR & other sources

DNA Clean & Concentrator™ (DCC™)

- ✓ Recovery of ultra-pure DNA that is free of salts and contaminants.
- ✓ Small (≥6 µl) elution volume.
- ✓ DNA is ideal for ligation, PCR, Next-Gen sequencing, etc.

Product	Size (Cat. No.)
DNA Clean & Concentrator™-5	50 Preps. (D4013) 200 Preps. (D4014)
ZR-96 DNA Clean & Concentrator™-5	2 x 96 Preps. (D4023) 4 x 96 Preps. (D4024)
Genomic DNA Clean & Concentrator™	25 Preps. (D4010) 100 Preps. (D4011)

High efficiency DNA recovery with the DCC™-5 compared to Supplier Q.

Boost DNA recoveries from agarose gels to >80%

Zymoclean™ Gel DNA Recovery

- ✓ Rapid (15 min.) recovery of ultra-pure DNA from agarose gels in ≥6 µl.
- ✓ Ultra-pure DNA ideal for DNA ligation, sequencing, etc.
- ✓ Format also available for large DNA >20 kb.


Product	Size (Cat. No.)
Zymoclean™ Gel DNA Recovery Kit	50 Preps. (D4001) 200 Preps. (D4002)
Zymoclean™ Large Fragment DNA Recovery Kit	25 Preps. (D4045) 100 Preps. (D4046)

DNA fragments recovered from an agarose gel using the Zymoclean™ Gel DNA Recovery Kit. Lanes: M: DNA Ladder; 1-5: individual ladder DNA fragments.

Recover transfection-quality plasmid DNA directly from culture

Zyppy™ Plasmid Prep Kits

- ✓ The fastest, simplest method available for purifying high quality plasmid DNA from *E. coli*.
- ✓ Pellet-Free™ procedure omits conventional cell-pelleting and resuspension steps.
- ✓ Transfection quality plasmid DNA directly from culture in under 15 minutes.

Product	Size (Cat. No.)
Zyppy™ Plasmid Miniprep Kit	50 Preps. (D4038)
	100 Preps. (D4019)
	400 Preps. (D4020)
	800 Preps. (D4037)

ZYMO RESEARCH CORP.

RNA PURIFICATION

Get RNA *directly* from TRIzol[®] without phase separation


What is Clean-Spin[™] Technology?
 The spin columns from Zymo Research have been designed to ensure complete elution with no binding/wash buffer carryover. The result is ultra-pure inhibitor-free DNA and RNA.

Direct-zol[™] RNA

- ✓ For purification of high-quality small and large RNA *directly* from TRIzol[®], TRI Reagent[®], or similar.
- ✓ Bypasses phase separation and precipitation procedures allowing for unbiased recovery of miRNA.

Product	Size (Cat. No.)
Direct-zol [™] RNA MiniPrep	50 Preps. (R2050)
	50 Preps. (R2051)*
	200 Preps. (R2052)
	200 Preps. (R2053)*

96-well and MagBead formats also available!
 DNase I included in all kits.
 * Supplied with TRI-Reagent[®]


High-quality small and large RNA are effectively recovered with the Direct-zol[™] kit. RNA is DNA-free.

Isolate DNA-free RNA from 1 to 10⁷ cells in minutes

Quick-RNA[™]

- ✓ Isolation of total, large, or small RNA – *You decide!*
- ✓ Ultra clean, high-quality RNA from a single cell to 10⁷ cells.
- ✓ DNA-free RNA ideal for any downstream application – *DNase I included.*

Product	Size (Cat. No.)
Quick-RNA [™] MicroPrep	50 Preps. (R1050)
	200 Preps. (R1051)
Quick-RNA [™] MiniPrep	50 Preps. (R1054)
	200 Preps. (R1055)
ZR-96 Quick-RNA [™]	2 x 96 Preps. (R1052)
	4 x 96 Preps. (R1053)

Isolate total, large, or small RNA with the Quick-RNA[™] kit. RNA is DNA-free using the Quick-RNA[™] kit.

Purify RNA from enzymatic and labeling reactions in 5 minutes

RNA Clean & Concentrator[™]

- ✓ Recover ultra-pure RNA in small ($\geq 6 \mu$ l) elution volumes.
- ✓ Compatible with TRIzol[®], phenol, chloroform, and RNase inhibitors (RNAlater[®]).
- ✓ RNA is ideal for RT-PCR, q-PCR, hybridization, arrays, RNA interference, etc.

Product	Size (Cat. No.)
RNA Clean & Concentrator [™] -5	50 Preps. (R1015)
	200 Preps. (R1016)
RNA Clean & Concentrator [™] -25	50 Preps. (R1017)
	100 Preps. (R1018)
ZR-96 RNA Clean & Concentrator [™]	2x96 well plates (R1080)
DNA-Free RNA Kit [™]	50 Preps. (R1013)
	200 Preps. (R1014)

The following are trademarks of other companies: pGEM[®], Promega Corp.; TRIzol[®] and TRI Reagent[®], Molecular Research Center, Inc.; Dri5[®] and Q1108[™], Life Technologies, Inc.

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Cycling Program for FRET Probes or Molecular Beacons

1 cycle of: 50°C for 10 min: Reverse transcription of RNA by M-MLV RT to generate first-strand cDNA.
1 cycle of: 95°C for 2 min: HotStart-IT Binding Protein and M-MLV RT inactivation.
35-45 cycles of: 95°C for 15 sec 50-60°C for 15-30 sec: Annealing temperature should be about 5°C below the T_m of the probe(s). Acquire real-time fluorescence data during this step.
72°C for 15-60 sec

Note: Due to the inherent specificity of TaqMan® probes, melt-curve analysis is not performed. Some probes such as FRET and Molecular Beacons or LUX™ primers may use melt-curve analyses. Consult the thermal-cycler manual and fluorescent probe/primer manufacturer recommendations for details.

7. If desired, confirm that specific PCR products have been generated by agarose gel electrophoresis. Amplicons may be detected on gels with ethidium bromide.

USB, the logo design, and HotStart-IT are registered trademarks of USB Corporation.

Tested User Friendly is a trademark of USB Corporation.

HotStart-IT Taq DNA Polymerase—Patent pending.

Taq DNA Polymerase—sold under licensing arrangements with Applied Biosystems. Purchase is accompanied by a limited license to use it in the Polymerase Chain Reaction (PCR) process in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Perkin-Elmer or as purchased, i.e., an authorized thermal cycler.

The Polymerase Chain Reaction (PCR) is covered by patents owned by Roche Molecular Systems and F. Hoffmann-La Roche Ltd.

TaqMan is a registered trademark of Roche Molecular Systems, Inc.

LUX is a trademark of Invitrogen.

ROX is a trademark of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

©USB Corporation 2008—All rights reserved.

USB Corporation

techsupport@usbweb.com

USB Europe GmbH

techsupport@usbweb.de

www.usbweb.com

P-75772B

Rev. 02/08

HotStart-IT® Probe One-Step qRT-PCR Master Mix Kit Tested User Friendly™ Product Number 75772

Brief Protocol

HotStart-IT® Probe One-Step qRT-PCR Master Mix Kit provides maximum convenience and optimal performance for real-time, quantitative analysis of RNA templates in a single reaction format. A separate tube of the passive reference dye, ROX, is included.

Protocol

This protocol applies to a single reaction where M-MLV RT, RNase Inhibitor, RNA template, primers, probe(s) and water are added to the HotStart-IT® Probe qPCR Master Mix. For multiple reactions, increase the volumes of the reaction components proportionally.

1. Thaw the master mix and other necessary frozen reagents at room temperature. Mix thoroughly, briefly spin to collect tube contents and then place on ice. The M-MLV RT and RNA samples should always be kept on ice.

2. Assemble reaction tubes or plates on ice.

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3. This table shows recommended component volumes. It is highly recommended to make a master mix for at least 10 reactions to reduce pipetting errors.

Component	Volume for 50 µl reaction	Volume for 20 µl reaction	Final Concentration
HotStart-IT [®] Probe qPCR Master Mix (2X)	25 µl	10 µl	1X (has 3.0mM MgCl ₂)
M-MLV RT	0.4 µl	0.16 µl	1X
RNase Inhibitor	0.4 µl	0.16 µl	1X
10µM Forward Primer	0.5-5.0 µl	0.2-2.0 µl	0.1-1.0µM*
10µM Reverse Primer	0.5-5.0 µl	0.2-2.0 µl	0.1-1.0µM*
10µM Fluorescent Probe(s)	0.5-2.5 µl	0.2-1.0 µl	0.1-0.5µM*
Template RNA	X µl	X µl	as needed, <1 µg†
ROX™ Passive Reference Dye (optional: for ABI and Stratagene instruments)	see included Dye protocols for details	Passive Reference	Reference
RNase-Free Water, DEPC-Treated	up to 50 µl	up to 20 µl	NA

*Optimal primer concentration is 0.2µM. In order to avoid primer-dimers and non-specific products, use ≤ 0.5µM. Because the reverse primer is also used during the initial reverse transcription step, it may be helpful to double the amount of reverse primer only. Start with a probe concentration of 0.1µM, although some optimization is generally required to determine the proper concentration.

†Total RNA may be used at 1 µg to 1 µg and poly(A)⁺ mRNA may be used at 100 fg to 100 ng per reaction.

NOTE: It is useful to test for contaminating genomic DNA in the RNA sample(s) by performing a control reaction which omits the M-MLV RT. Since there is no reverse transcriptase in this reaction, any observed product must have been generated from a DNA source. Also, to confirm the absence of general DNA or RNA contamination in other reagents, perform a No-Template-Control reaction.

4. *Optional:* If optimizing Mg²⁺ concentration, add 2.0 µl of 25mM MgCl₂ per 50 µl reaction for each additional 1mM Mg²⁺ required. Subtract this volume from the amount of water needed.

5. Cap tubes or seal plates with optically clear caps or film. Mix tubes or plates by gentle vortexing and then spin to collect contents without bubbles (e.g. 2-5 min at 1000-2000 x g).

6. The following tables show recommended cycling conditions:

Standard Cycling Program for TaqMan[®] Probes

1 cycle of: 50°C for 10 min: Reverse transcription of RNA by M-MLV RT to generate first-strand cDNA.
1 cycle of: 95°C for 2 min: HotStart-IT Binding Protein and M-MLV RT inactivation.
35-45 cycles of: 95°C for 15 sec 60°C for 30-60 sec: Acquire real-time fluorescence data during this step.

Fast Cycling Program (e.g. ABI 7500 in Fast Mode)

1 cycle of: 50°C for 5 min: Reverse transcription of RNA by M-MLV RT to generate first-strand cDNA.
1 cycle of: 95°C for 2 min: HotStart-IT Binding Protein and M-MLV RT inactivation.
35-45 cycles of: 95°C for 3-5 sec 60°C for 15-30 sec: Acquire real-time fluorescence data during this step.

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